



Dear customer

23 June 2016

ULVAC Components News

List Price Revision of Ionization Vacuum Gauge Sensor Head WIT,WIB and WIN

We regret to inform that we will revise price of following sensor heads for glass type and nude type vacuum gauges. We had kept current list price for more than 40 years in spite of frequent cost increase of manufacture reduction. However, there is no way anymore to sustain by only our management effort. We thank you for your understanding.

We appreciate your business and look forward to serving you in the future.

- Detail -

(1) Effective date of new price

From the purchase on January 4th, 2017.

(2) Corresponding Model

Sensor head for ionization vacuum gauge (WIT,WIB,WIN)

Price revision will be applied only to the sensor head, the price of Ionization vacuum be considered as stationary.

(GI-D7 is considered in stationary price)

(3) New price list

| Model | Specifications | New List price (Yen) | Old List price (Yen) |
|-------------------------|--------------------|----------------------|----------------------|
| Sensor head WIT-G1 1pc | Φ 18/Φ 15 | 15,000 | 10,500 |
| Sensor head WIT-G8 1pc | Φ 18/Φ 15 | 48,000 | 40,200 |
| Sensor head WIB-G5 1pc | Φ 18/Φ 15 | 21,000 | 15,000 |
| Sensor head WIB-N8 1pc | Nude Type | 68,000 | 60,000 |
| | | | |
| Sensor head WIN-N1 1pc | Without flange | 98,000 | 85,000 |
| Sensor head WIN-N2 1pc | Pin type connector | 98,000 | 87,000 |
| Sensor head WIN-N3 1pc | BNC Type connector | 125,000 | 108,000 |
| Filament WIN-N Type 1pc | Tungsten | 16,000 | 13,000 |

Sincerely

TAKUYA KOBAYASHI

Product Manager Group Dept.2 Senior Manager

ULVAC Inc. Components Div.

End of Document