



Dear customers,

October 2017

ULVAC Components News

End of support and maintenance

Please be informed that we stopped supporting following obsolete product due to termination of support and maintenance term.

Please consider following successor model in the future.

– Detail –

The Sales support of Turbo Molecular Pump for HELIOT 300 series will end at the end of December 2018.

We continue to maintain HELIOT 300 series as long as possible. However, there will be unavailable to supply due to the end production.

Furthermore, in order to prevent troubles caused by the inability to repair HELIOT 300 we recommend the upgrade of HELIOT 300 series to HELIOT 900 series.

1. Object model/Maintenance term
Leak detector HELIOT300series / Turbo Molecule Pump
2. Succeeding model
Leak detector HELIOT900series

Sincerely

TAKUYA KOBAYASHI
Product Manager Group Dept.2 Senior Manager
ULVAC Inc. Components Div.

End of Document